



WandelWol®

The solution for blisters and foot discomfort



100% natural product



Suitable for all sports



Easy to apply



Supports regenerative abilities of the skin.

WandelWol anti-pressure wool

The anti-pressure product by WandelWol is made from 100% sheep wool; without chemical additions or modifications. The material is soft, contains lanolin (wool wax) and makes you feel like you are walking on a cloud.

Because this product is made from natural fibres, it is biodegradable. WandelWol products absorb perspiration and allows for ventilation. It also helps keep your feet warm in winter and helps cool them in summer. The anti-pressure wool softens blisters and prevents new ones from forming. It also takes away pressure from the feet and helps preserve toenails and feet.

For who?

WandelWol anti-pressure products are suitable for hikers, athletes, shoppers and people that have little padding under their feet, cold feet, painful calluses in between or around the toes and with a burning sensation under the ball of the foot. Prevents the forming of trauma nails during exercise or long hikes. For instructions check the back of this sheet.



For more information see www.wandelwol.nl or send an email to info@wandelwol.nl

Instructions

1. When you sense a painful spot, stop what you are doing and apply the product. Do not wait for a blister to form.
2. Apply the WandelWol anti-pressure wool directly to the painful spot and put on your socks and shoes like normal. It will give relief immediately. Depending on the spot, you might not need to take off your sock completely to apply.
3. If the problem occurs around your toes, remove the sock and wrap the anti-pressure wool around the painful area and put your sock on over the wool.
4. WandelWol products integrate into the sock over the day, this is normal. The removal of the anti-pressure wool is very easy. Change it out if necessary.
5. WandelWol product retain the developed fit to the foot. The anti-pressure wool does not form lumps or roll up.
6. After using it a few times you will find the right amount of product for you.
7. WandelWol anti-pressure wool is sold in bags of 10, 20 and 40 grams with a set of instructions inside.
8. WandelWol anti-pressure 10 grams is suitable for locations in between the toes, 20 grams is for larger locations on the foot or to relieve larger toes and a bag of 40 grams is suitable for when you hike or exercise a lot.
9. Suitable for all sports and best for hiking, running, playing football etc.
10. Use for relieving (large) toes to prevent trauma nails.
11. Use for preventing and relieving blisters, pain relief from calluses (it is advised to have these removed by a licenced pedicurist), burning and tired feet.
12. WandelWol products support the regenerative abilities of the skin.

WandelWol anti-pressure products should not be used on damaged skin or open wounds - it is not medication.

Do not use with peripheral neuropathy (numbness) - patients should ask a doctor before using the product.

Only use WandelWol products temporarily - it is not a long-term solution! Visit your pedicurist or podiatrist regularly when suffering from chronic pain.

BE AWARE: If allergic to lanolin (wool wax), you should not use WandelWol anti-pressure wool based on sheep wool. It is advised to use WandelWol anti-pressure wool based on alpaca fiber as this does not contain lanolin.



0-3